


CRYOCARE SL™

The Cryocare SL™ system is a streamlined version of the Cryocare CS®, designed specifically for the percutaneous ablation market. Cryocare SL delivers the same high level of performance and ergonomic functionality as the CS system, however, with a simpler more streamlined operational and software interface.



19" Color Monitor (Articulation 315°)

3 way locking display height. Max height 24"

Adjustable viewing lock

Audio Output

Braking wheel lock

Removable Remote Keypad

Metal lid for storage protection

One piece handle

Rewritable CD Drive

Space saving, ergonomic base

Directional wheel lock

Cryocare SL Specifications

System Size (L x W x H)	Approximately 30" x 24" x 40"
Weight	<300 lbs.
Electrical	120/230 VAC 50/60 Hz 4 Amps RMS
Fuse	T6 3A 250V

 **endocare™**
extending life everyday™
800-418-4677