


CRYOCARE[™] CS

with AutoFreeze

Cryocare CS is a fully integrated planning, placement and treatment system designed to simplify cryotherapy and allow unlimited flexibility. Cryocare CS is designed to meet specific physician needs through innovative, user-friendly technology that introduces a new generation of cryotherapy.



AutoFreeze[™] Software

20" Color Monitor
(Articulation 315°)

3 way locking display
height. Max height 24"

Adjustable
viewing lock

2 way removable
extending probe holder

Color backlit
for low light use

Ultrasound with
Color Doppler

90 degree articulating
& foldable control panel

Onboard Printer

Audio Output

Foot pedal to
change U/S view

Urethral Warmer

Removable
Remote Keypad

Metal lid for
storage protection

BNC video output

One piece handle

Integrated Ultrasound

Rewritable CD Drive

Space saving,
ergonomic base

Braking wheel lock

Directional
wheel lock

Cryocare CS Specifications

System Size (L x W x H)	Approximately 30" x 24" x 40"
Weight	<300 lbs.
Electrical	120/230 VAC 50/60 Hz 4 Amps RMS
Fuse	T6 3A 250V

Integrated Ultrasound Specifications

Ultrasound probe frequencies	5-10 MHz (10L5)
Frame rate	124 Hz max (based on incoming data rate)
Ultrasound lines/frame	128

 **endocare[™]**
extending life everyday[™]
800-418-4677